

PKIG Antibody



PACO11270

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PKIG

Uniprot

Q9Y2B9

Synonyms:

protein kinase (cAMP-dependent, catalytic) inhibitor gamma; PKIG; MGC126458; MGC126459 ;

Immunogen:

Human PKIG.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A