

PKLR Antibody



PACO11271

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PKLR

Uniprot

P30613

Synonyms:

pyruvate kinase, liver and RBC; PKLR; PK1; PKL; PKR; PKRL; RPK ;

Immunogen:

Human PKLR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A