

PKNOX2 Antibody



PACO11276

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PKNOX2

Uniprot

Q96KN3

Synonyms:

PBX/knotted 1 homeobox 2; PKNOX2; FLJ13074; PREP2 ;

Immunogen:

Human PKNOX2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A