

PLCH1 Antibody



PACO11289

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PLCH1

Uniprot

Q4KWH8

Synonyms:

phospholipase C, eta 1; PLCH1; DKFZp434C1372; MGC117152; PLCL3; PLCeta1 ;

Immunogen:

Human PLCH1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A