

# BLOC1S6 Antibody



PACO11292

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BLOC1S6

**Uniprot**

Q9UL45

**Synonyms:**

pallidin homolog (mouse); PLDN; PA; PALLID ;

**Immunogen:**

Human PLDN.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A