

PLEKHA1 Antibody



PACO11294

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PLEKHA1

Uniprot

Q9HB21

Synonyms:

pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 1; PLEKHA1; TAPP1 ;

Immunogen:

Human PLEKHA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A