

PLEKHA3 Antibody



PACO11295

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PLEKHA3

Uniprot

Q9HB20

Synonyms:

pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 3; PLEKHA3; FAPP1; FLJ20067 ;

Immunogen:

Human PLEKHA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A