

PLEKHF1 Antibody



PACO11299

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PLEKHF1

Uniprot

Q96S99

Synonyms:

pleckstrin homology domain containing, family F (with FYVE domain) member 1; PLEKHF1; APPD; MGC4090; PHAFIN1; ZFYVE15 ;

Immunogen:

Human PLEKHF1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A