

PLEKHG5 Antibody



PACO11300

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PLEKHG5

Uniprot

O94827

Synonyms:

pleckstrin homology domain containing, family G (with RhoGef domain) member 5; PLEKHG5; DSMA4; GEF720; KIAA0720; RP4-650H14.3 ;

Immunogen:

Human PLEKHG5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A