

PLP2 Antibody



PACO11307

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PLP2

Uniprot

Q04941

Synonyms:

proteolipid protein 2 (colonic epithelium-enriched); PLP2; A4; A4-LSB; MGC126187 ;

Immunogen:

Human PLP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A