

PLXDC2 Antibody



PACO11313

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PLXDC2

Uniprot

Q6UX71

Synonyms:

plexin domain containing 2; PLXDC2; FLJ14623; TEM7R ;

Immunogen:

Human PLXDC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A