

PM20D2 Antibody



PACO11314

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PM20D2

Uniprot

Q8IYS1

Synonyms:

peptidase M20 domain containing 2; PM20D2; ACY1L2; bA63L7.3 ;

Immunogen:

Human PM20D2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A