

PACO11324

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PNCK

Uniprot

Q6P2M8

Synonyms:

pregnancy up-regulated non-ubiquitously expressed CaM kinase; PNCK; BSTK3; CaMK1b; FLJ50403; FLJ50549; FLJ56451; FLJ59811; MGC45419 ;

Immunogen:

Human PNCK.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A