

PNMA1 Antibody



PACO11327

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PNMA1

Uniprot

Q8ND90

Synonyms:

paraneoplastic antigen MA1; PNMA1; MA1 ;

Immunogen:

Human PNMA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A