

PNMT Antibody



PACO11330

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PNMT

Uniprot

P11086

Synonyms:

phenylethanolamine N-methyltransferase; PNMT; MGC34570; PENT; PNMTase ;

Immunogen:

Human PNMT.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A