

PNPLA6 Antibody



PACO11333

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PNPLA6

Uniprot

Q8IY17

Synonyms:

patatin-like phospholipase domain containing 6; PNPLA6; NTE; NTEMND; SPG39; sws ;

Immunogen:

Human PNPLA6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A