

# PNRC1 Antibody



PACO11334

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

PNRC1

**Uniprot**

Q12796

**Synonyms:**

proline-rich nuclear receptor coactivator 1; PNRC1; B4-2; PNAS-145; PROL2; PRR2; RP11-63L7.5 ;

**Immunogen:**

Human PNRC1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A