

PODXL2 Antibody



PACO11336

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PODXL2

Uniprot

Q9NZ53

Synonyms:

podocalyxin-like 2; PODXL2; PODLX2 ; Endoglycan;

Immunogen:

Human PODXL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A