

# POLA2 Antibody



PACO11340

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POLA2

**Uniprot**

Q14181

**Synonyms:**

polymerase (DNA directed), alpha 2 (70kD subunit); POLA2; FLJ21662; FLJ37250 ;

**Immunogen:**

Human POLA2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A