

POLR1A Antibody



PACO11351

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

POLR1A

Uniprot

O95602

Synonyms:

polymerase (RNA) I polypeptide A, 194kDa; POLR1A; DKFZp586M0122; FLJ21915; MGC87965; RPA1; RPO1-4 ; RPA194

Immunogen:

Human POLR1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A