

POLR1C Antibody



PACO11352

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

POLR1C

Uniprot

O15160

Synonyms:

polymerase (RNA) I polypeptide C, 30kDa; POLR1C; RPA39; RPA40; RPA5; RPAC1 ;

Immunogen:

Human POLR1C.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A