

POLR2C Antibody



PACO11355

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

POLR2C

Uniprot

P19387

Synonyms:

polymerase (RNA) II (DNA directed) polypeptide C, 33kDa; POLR2C; RPB3; RPB31; hRPB33; hsRPB3 ;

Immunogen:

Human POLR2C.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A