

POLR2D Antibody



PACO11356

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POLR2D

Uniprot

O15514

Synonyms:

polymerase (RNA) II (DNA directed) polypeptide D; POLR2D; HSRBP4; HSRPB4; RBP4 ;

Immunogen:

Human POLR2D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A