

POLR2E Antibody



PACO11357

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

POLR2E

Uniprot

P19388

Synonyms:

polymerase (RNA) II (DNA directed) polypeptide E, 25kDa; POLR2E; RPABC1; RPB5; XAP4; hRPB25; hsRPB5 ;

Immunogen:

Human POLR2E.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A