

POLR2F Antibody



PACO11358

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POLR2F

Uniprot

P61218

Synonyms:

polymerase (RNA) II (DNA directed) polypeptide F; POLR2F; HRBP14.4; POLRF; RPABC2; RPB14.4; RPB6 ;

Immunogen:

Human POLR2F.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A