

POLR3A Antibody



PACO11362

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POLR3A

Uniprot

O14802

Synonyms:

polymerase (RNA) III (DNA directed) polypeptide A, 155kDa; POLR3A; RPC1; RPC155; hRPC155 ;

Immunogen:

Human POLR3A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A