

# POLR3B Antibody



PACO11363

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

POLR3B

**Uniprot**

Q9NW08

**Synonyms:**

polymerase (RNA) III (DNA directed) polypeptide B; POLR3B; FLJ10388; RPC2 ;

**Immunogen:**

Human POLR3B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A