

POP1 Antibody



PACO11371

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POP1

Uniprot

Q99575

Synonyms:

processing of precursor 1, ribonuclease P/MRP subunit (S. cerevisiae); POP1; KIAA0061; MGC17365 ;

Immunogen:

Human POP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A