

# PPM1L Antibody



PACO11412

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PPM1L

**Uniprot**

Q5SGD2

**Synonyms:**

protein phosphatase 1 (formerly 2C)-like; PPM1L; MGC132545; MGC132547; PP2C-epsilon; PP2CE; PPM1-LIKE ;

**Immunogen:**

Human PPM1L

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A