

PPP1CB Antibody



PACO11416

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PPP1CB

Uniprot

P62140

Synonyms:

protein phosphatase 1, catalytic subunit, beta isoform; PPP1CB; MGC3672; PP-1B; PPP1CD ;

Immunogen:

Human PPP1CB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A