

# PPP1CC Antibody



PACO11417

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

PPP1CC

**Uniprot**

P36873

**Synonyms:**

protein phosphatase 1, catalytic subunit, gamma isoform; PPP1CC; PPP1G ;

**Immunogen:**

Human PPP1CC.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A