

PPP1R12B Antibody



PACO11418

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PPP1R12B

Uniprot

O60237

Synonyms:

protein phosphatase 1, regulatory (inhibitor) subunit 12B; PPP1R12B; MGC131980; MGC87886; MYPT2 ;

Immunogen:

Human PPP1R12B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A