

PRICKLE3 Antibody



PACO11461

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PRICKLE3

Uniprot

O43900

Synonyms:

prickle homolog 3 (Drosophila); PRICKLE3; LMO6 ;

Immunogen:

Human PRICKLE3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A