

PRPF19 Antibody



PACO11495

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PRPF19

Uniprot

Q9UMS4

Synonyms:

PRP19/PSO4 pre-mRNA processing factor 19 homolog (S. cerevisiae); PRPF19; NMP200; PRP19; PSO4; SNEV; UBOX4; hPSO4 ;

Immunogen:

Human PRPF19.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A