

PRPF40A Antibody



PACO11498

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PRPF40A

Uniprot

O75400

Synonyms:

PRP40 pre-mRNA processing factor 40 homolog A (S. cerevisiae); PRPF40A; FBP-11; FBP11; FLAF1; FLJ20585; FBNP3; HIP10; HYPA; NY-REN-6 ;

Immunogen:

Human PRPF40A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A