

# PRUNE2 Antibody



PACO11513

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRUNE2

**Uniprot**

Q8WUY3

**Synonyms:**

prune homolog 2 (Drosophila); PRUNE2; A214N16.3; BMCC1; BNIPXL; C9orf65; DKFZp762K117; FLJ50060; FLJ54876; FLJ59118; KIAA0367; RP11-58J3.2; bA214N16.3 ;

**Immunogen:**

Human PRUNE2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A