

PSG11 Antibody



PACO11520

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PSG11

Uniprot

Q9UQ72

Synonyms:

pregnancy specific beta-1-glycoprotein 11; PSG11; MGC22484; PSG13; PSG14 ;

Immunogen:

Human PSG11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A