

PSG4 Antibody



PACO11521

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PSG4

Uniprot

Q00888

Synonyms:

pregnancy specific beta-1-glycoprotein 4; PSG4; PSG9 ; CGM4, PSG9

Immunogen:

Human PSG4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A