

# PSMB10 Antibody



PACO11530

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PSMB10

**Uniprot**

P40306

**Synonyms:**

proteasome (prosome, macropain) subunit, beta type, 10; PSMB10; LMP10; MECL1; MGC1665; beta2i ;

**Immunogen:**

Human PSMB10.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A