

PSMB11 Antibody



PACO11531

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PSMB11

Uniprot

A5LHX3

Synonyms:

proteasome (prosome, macropain) subunit, beta type, 11; PSMB11; FLJ16369; beta5t ;

Immunogen:

Human PSMB11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A