

PSMD10 Antibody



PACO11545

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PSMD10

Uniprot

O75832

Synonyms:

proteasome (prosome, macropain) 26S subunit, non-ATPase, 10; PSMD10; dJ889N15.2; p28 ;

Immunogen:

Human PSMD10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A