

PTGR2 Antibody



PACO11575

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PTGR2

Uniprot

Q8N8N7

Synonyms:

prostaglandin reductase 2; PTGR2; DKFZp686P10120; FLJ39091; FLJ99229; PGR2; ZADH1 ;

Immunogen:

Human PTGR2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A