

# PTMA Antibody



PACO11581

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PTMA

**Uniprot**

P06454

**Synonyms:**

prothymosin, alpha; PTMA; MGC104802; TMSA ;

**Immunogen:**

Human PTMA.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A