

# PTPRCAP Antibody



PACO11595

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PTPRCAP

**Uniprot**

Q14761

**Synonyms:**

protein tyrosine phosphatase, receptor type, C-associated protein; PTPRCAP; CD45-AP; LPAP; MGC138602; MGC138603 ;

**Immunogen:**

Human PTPRCAP.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A