

PTTG1IP Antibody



PACO11607

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PTTG1IP

Uniprot

P53801

Synonyms:

pituitary tumor-transforming 1 interacting protein; PTTG1IP; C21orf1; C21orf3; PBF ;

Immunogen:

Human PTTG1IP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A