

PUF60 Antibody



PACO11609

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PUF60

Uniprot

Q9UHX1

Synonyms:

poly-U binding splicing factor 60KDa; PUF60; FIR; FLJ31379; RoBPI; SIAHBP1 ;

Immunogen:

Human PUF60.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A