

PUS3 Antibody



PACO11613

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PUS3

Uniprot

Q9BZE2

Synonyms:

pseudouridylate synthase 3; PUS3;2610020J05Rik; FKSG32; FLJ23638 ;

Immunogen:

Human PUS3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A