

RAB11A Antibody



PACO11630

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RAB11A

Uniprot

P62491

Synonyms:

RAB11A, member RAS oncogene family; RAB11A; MGC1490; YL8 ;

Immunogen:

Human RAB11A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A