

RAB12 Antibody



PACO11633

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RAB12

Uniprot

Q6IQ22

Synonyms:

RAB12, member RAS oncogene family; RAB12; FLJ45927; MGC104724 ;

Immunogen:

Human RAB12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A