

# RAB20 Antibody



PACO11638

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RAB20

**Uniprot**

Q9NX57

**Synonyms:**

RAB20, member RAS oncogene family; RAB20; FLJ20429 ;

**Immunogen:**

Human RAB20.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A